


EGIS Flower Stent

Upgraded Special Features



- Innovative design with 5~6 grooves, facilitating bile drainage



- Fully covered stent with e-PTFE, designed for use in malignant & benign strictures



- Sufficient radial force to push out the stricture



- Low axial force, ensuring excellent flexibility & conformability



- For easy removal, a triangle-shaped lasso is attached



- To obtain anti-migration function, flared ends and grooves are made through out the entire body

Specifications

Delivery System 8Fr, 50/180cm

Order NO. First two** : Diameter(mm) / Last two** : Length(cm)

Order No.			Size	
Type	OTW	TTS	Diameter	Length
Bare	BFBP****0805	BFBP****0818	8, 10	4, 5, 6, 7, 8, 9, 10, 12
Partial	BFPP****0805	BFPP****0818		
Fully Covered	BFCP****0805	BFCP****0818		

