



SEAL THORACIC STENT GRAFT

Inward Bending
Excellent Proximal Radial Force
Complete Sealing Zone
Low Graft Permeability
Excellent Flexibility & Conformability
4, 6, 8, 10mm Tapered




CASE



Case1

Pre CT 3D image

Post CT 3D image




Case2


Pre CT 3D image

Post CT 3D image


Graft




SEAL Thoracic Stent Graft



Competitor A



SEAL Thoracic Stent Graft
(Magnified appearance)



Comparison Chart

Diameter (mm)	Vasnova Delivery system O.D	A Company Delivery system O.D	B Company Delivery system O.D
	24	18Fr (5.92mm)	18Fr
26			
28			
30			
32			
34	21Fr (7.0mm)	21Fr	24Fr
36			
38			
40	21Fr (7.0mm)	23Fr	25Fr
42			
44			
46			

Specification

Type	Order No.	Size			
		Diameter (mm)	Length (mm)	Delivery System	
				Diameter OD(ID)Fr	Usable Length(cm)
Thoracic Limb Stent	Limb-**-** Limb-**-***	22	60 ~ 100 110 ~ 200	15 (13)	90
		24 ~ 36		18 (16)	
		38 ~ 46		21 (19)	

Order No. 1st **: Diameter(mm) / Last **: Length(mm)

Type	Order No.	Diameter (mm)		Length (mm)	Delivery System	
		Proximal Diameter (mm)	Distal Diameter (mm)		Diameter OD(ID)Fr	Usable Length(cm)
		Thoracic Tapered Limb Stent (4, 6, 8, 10 mm)	Limb-**-*** Limb-**-****	24	20	80, 110, 130, 150, 180, 200
26	20,22					
28	20,22,24					
30	20,22,24,26					
32	22,24,26,28					
34	24,26,28,30					
36	26,28,30,32			80, 110, 130, 150, 180, 200	21 (19)	90
38	28,30,32,34					
40	30,32,34,36					
42	32,34,36,38					
44	34,36,38,40					
46	36,38,40,42					

Order No. 1st **: Proximal Diameter(mm) / 2nd **: Distal Diameter(mm) / 3rd **: Length(mm)