


EGIS Colorectal Stent

Rectal Stent for large bowel obstruction



- ▣ Sufficient radial force
- ▣ Outstanding flexibility and conformability without kinking
- ▣ Excellent visibility under fluoroscopy

Specification

Delivery System 10Fr, 80/150/220cm | 12Fr, 80/150cm

Ordering NO. First two** : Diameter(mm) / Last two** : Length(cm)

TYPE	Stent		Delivery System	
	Diameter	Length	Profile(Fr)	Length(cm)
Single Bare	18, 20, 22, 24	6, 8, 10, 12, 14	10	80/150/220
	26, 28, 30		12	80/150
Single Covered	18, 20, 22, 24	6, 8, 10, 12	10	80/150/220
Double Bare	26, 28, 30		12	80/150
Double Covered				

- ▣ Clinical success rate of EGIS covered and uncovered colorectal stents was 97%, re-obstruction was rare in both groups.
- ▣ New design Double Bare or Double Covered Stents for Malignant Colorectal Obstruction:
Comparative, Prospective, Randomized, Single Centre Trial
A Vodoleev, V Duvanskiy, D Kryazhev Endoscopy 2018; 50(04): S33



S&G BIOTECH INC.
www.sng-bio.com


82, Bugok-ro, Pogok-eup, Cheoin-gu, Yongin-si, Gyeonggi-do, Korea (17023)
e-mail. sng@sng-bio.com Tel.+82 31.748.6625 Fax.+82 31.748.6620




CAD-03-01(Rev.0)

EGIS Pyloric Stent

Duodenal Stent for gastric outlet obstruction



- Sufficient radial force
- Outstanding flexibility and conformability without kinking
- Excellent visibility under fluoroscopy



■ EGIS Pyloric Double Covered Stent



■ EGIS Pyloric Double Bare stent

Specification

Delivery System 10Fr, 120/180cm

Order NO. First two** : Diameter(mm) / Last two** : Length(cm)

Type	OrderNO.		Size		Delivery system (Fr)
	OTW	TTS	Diameter	Length	
Single Bare	SEBP****	SEBE****	18, 20 22, 24	6, 8, 10, 12, 14	10
Single Covered	SECP****	SECE****		6, 8, 10, 12	
Double Bare	DEBP****	DEBE****		6/5, 8/7, 10/9, 12/11	
Double Covered	DECP****	DECE****			



S&G BIOTECH INC.
www.sng-bio.com

82, Bugok-ro, Pogok-eup, Cheoin-gu, Yongin-si, Gyeonggi-do, Korea (17023)
e-mail. sng@sng-bio.com Tel.+82 31.748.6625 Fax.+82 31.748.6620



CAD-04-01(Rev.0)