


EGIS Biliary Stent Braided 6



01 Single Bare stent




- Excellent flexibility & conformability
- Small cell size prevents tumor in-growth.
- Low axial force

Delivery system 6/8Fr, 50/180cm
Order No. First two** : Diameter(mm) / Last two** : Length(cm)

Order No.		Size		Delivery system
OTW	TTS	Diameter(mm)	Length(cm)	Diameter(Fr)
SBP****06F	SBE****06F	6, 7, 8, 10	3, 4, 5, 6, 7, 8, 9, 10, 12	6
SBP****0	SBE****0	8, 10, 12	4, 5, 6, 7, 8, 9, 10, 12	8
KSB****05	KSB****18			

02 Single Covered stent




- Fully covered with silicone membrane.
- For easy removal, lasso is attached.

Delivery system 6/8Fr, 50/180cm
Order No. First two** : Diameter(mm) / Last two** : Length(cm)

Order No.		Size		Delivery system
OTW	TTS	Diameter(mm)	Length(cm)	Diameter(Fr)
SCP****0	SCE****0	6	3, 4, 5, 6, 7, 8, 9, 10	6
KSC****05	KSC****18	8, 10	4, 5, 6, 7, 8, 9, 10, 12	8

03 Double Bare stent




- Excellent flexibility & conformability
- Crossly overlapped two stents provide smaller cell size reducing the rate of tumor in-growth.
- Sufficient radial force

Delivery system 6/8Fr, 50/180cm
Order No. First two** : Diameter(mm) / Last two** : Length(cm)

Order No.		Size		Delivery system
OTW	TTS	Diameter(mm)	Length(cm)	Diameter(Fr)
DBP****0	DBE****0	7	4, 5, 6, 7, 8, 9, 10	6
KDB****05	KDB****18	8, 10	4, 5, 6, 7, 8, 9, 10, 12	8

04 SCF cone shape stent



- Tapered flower design
- Full covered silicone membrane
- For Hiliar lesion

Delivery system 6Fr, 180cm
Order No. two** : Length(cm)

Order No.	Size		Lasso
	Diameter(mm)	Length(cm)	Length(cm)
SECE07**0FC	7/6	4, 5, 6, 7, 8, 9, 10, 12, 14	10, 9, 8, 7, 6, 4, 2

