

EGIS Biliary Partial DC 20 Stent





Superior silicone sleeve placed between two layers of stents

Excellent flexibility & conformability

Smaller cell size reduces the rate of tumor-ingrowth

Prevents side branch obstruction



Specification

Delivery System 8Fr, 500/1800mm
 Order No. First two** : Diameter(mm) / Last two** : Length(mm)

Type	Order information			
	Order No.		Size	
	OTW	TTS	Diameter	Length
Partial DC 20 Stent	DPP****0	DPE****0	8, 10	60, 70, 80, 90, 100, 120



S&G BIOTECH INC.
 www.sng-bio.com

82, Bugok-ro, Pogok-eup, Cheoin-gu, Yongin-si, Gyeonggi-do, Korea (17023)
 e-mail, sng@sng-bio.com Tel,+82 31,748,6625 Fax,+82 31,748,6620



CAD-01-02(Rev.0)