




# RadioMate II

Microcatheter

Super **Selective** and  
High **Flow** are in your choice  
Being an Innovative Leading Edge




VASNOVA


## Features & Benefits

- ✔ **Hydrophilic Coating**  
 Gives better trackability through tortuous anatomy
- ✔ **Stainless steel dual coil braiding**  
 Increases kink resistance, making it possible to approach to all peripheral and visceral vessels
- ✔ **3-layer structure**  
 Protects inner lumen integrity for smooth delivery of all types of embolic materials

### 3-Layer shaft structure



## ORDER INFORMATION



Order No.	O.D. of distal shaft(Fr)	O.D. of proximal shaft(Fr)	Inner diameter(mm)	Max. GW O.D.(inch)	Distal Marker	Usable Length(cm)
JG-HMC-***	1.6	2.6	0.38	0.014	0.6mm length /0.6mm from distal end	110, 120, 130, 140, 150
JG-KMC-***	1.8	2.7	0.41	0.016	0.6mm length /0.6mm from distal end	
JG-JMC-***	2.1	2.7	0.53	0.016	0.6mm length /0.6mm from distal end	
JG-YMC-***	2.3	2.7	0.54	0.016	0.6mm length /0.6mm from distal end	
JG-WMC-***	2.7	2.7	0.69	0.018	0.6mm length /0.6mm from distal end	110, 120, 130

(Order No. \*\*\* : Length (mm) )

Hydrophilic coating length : 60cm