

EGIS Biliary Wide Stent



stent BARE



2cm


Longer bare
Anti-migration

Features & Benefits

- ▣ Reduces infection significantly
- ▣ Excellent dependability, confirmability with a double layer structure.
- ▣ High biocompatible e-PTFE membrane

stent DIAMETER

12mm



44% increased inner volume

Specification

Delivery system : 8Fr, 50/180cm

Order NO. First two** : Diameter(mm) / Last two** : Length(cm)

Order No.		Size	
OTW	TTS	Diameter(mm)	Length(cm)
DCP****OE	DCE****OE	12	4, 5, 6, 7, 8, 9, 10