


SEAL Sheath Dilator

Accessory for SEAL Stent graft



Guiding Sheath For stent Graft

1. Minimize blood loss with double locking hemostatic valve system
2. Use a large size of the sheath in the SEAL stent graft placement
3. The excellent flexibility with sheath braiding treatment

Ordering information

Ordering No.	Sheath		Dilator	
	Diameter(Fr)	Usable Length(mm)	Diameter(mm)	Length(mm)
N01092140300	14	300	4.76	416
N01092160300	16	300	5.44	414
N01092180300	18	300	6.07	414
N01092200300	22	300	7.44	424
N01092240350	24	350	8.08	484

