

Pliant Guide wire



Hydrophilic coating

Core-to-tip design

Nitinol alloy core

Polyurethane jacket with tungsten

Features and Benefits

Hydrophilic coating

- Unmatched lubricity consistently gets you where you need to go

Core-to-tip design

-Provides optimal torque control for excellent navigation

Polyurethane jacket with tungsten

-Provide enhanced radiopacity of entire guide wire

Nitinol alloy core


-Resists kinking for easier, faster navigation through even the most tortuous vessels

Multi-Purpose Guide wire

-Extra-Regular type Guide wire used for all procedures of the Intervention

Distal Radiopaque Point

-Enhancement of visibility on the fluoroscopy



Ordering Information

Tip type	Usable length(cm)	Ordering No.	Diameter(inch)	Flexible tip length(cm)
Angled	80	PGA-35-80	0.035	3
	100	PGA-35-100		
	120	PGA-35-120		
	150	PGA-35-150		
	180	PGA-35-180		
	260	PGA-35-260		
Straight	80	PGS-35-80	0.035	3
	100	PGS-35-100		
	120	PGS-35-120		
	150	PGS-35-150		
	180	PGS-35-180		
	260	PGS-35-260		

