



EGIS Biliary Single Covered F Stent



Features & Benefits

- ▣ Creative design for smooth flow of bile through the tubes resistant to compression due to 5 or 6 groove shape
- ▣ Distal CBD optimization through expansion forces improves pancreatic drainage
- ▣ Anti-migration due to Silicone membrane



Specification

Delivery system : 8Fr, 50/180cm

Order NO. First two** : Diameter(mm) / Last two** : Length(cm), Tolerance:±5%

Order No.		Size	
OTW	TTS	Diameter	Length
SFP****0	SFE****0	8,10	4, 5, 6, 7, 8, 9, 10, 12

* Removal Lasso Length : 2CM



S&G BIOTECH INC.
www.sng-bio.com

82, Bugok-ro, Pogok-eup, Cheoin-gu, Yongin-si, Gyeonggi-do, Korea (17023)
e-mail. sng@sng-bio.com Tel.+82 31.748.6625 Fax,+82 31.748.6620



CAD-01-12(Rev.0)