


Bifurcation Support
Kink resistance


SEAL NOVUS STENT GRAFT

Complete Sealing Zone
Low Graft Permeability
Juxtarenal Fixation
6 Barbs for anti-migration


SEAL NOVUS Body Stent Graft



Inward Bending &
6 Barbs for anti-migration



Running Suture



Endo-Skeleton &
Diamond Cell Structure




Bifurcation Support
Cell Structure

CUFF Stent Graft



NOVUS LIMB Stent Graft



NOVUS Flared LIMB Stent Graft



Specification

Type	Stent		Delivery System	
	Diameter(mm)	Length(mm)	Out Diameter(Fr)	Length(mm)
SEAL NOVUS Body Stent Graft	22	50	15	550
	24,26,28,30		18	
	32,34,36		21	
NOVUS Limb Stent Graft	12	40,60,80,100,120	15	
NOVUS Flared Limb Stent Graft	12/14,16,18,20,22,24	60,80,100,120	15	
CUFF Stent Graft	20,22	40	15	
	24,26,28,30, 32,34,36,38		18	