






# EGIS Biliary Hilar FLEX



- The Hilar FLEX has a structure that **allows cell expansion only in the central part** to enable Stent-In-Stent procedures.
- To make it easier for the delivery system to enter the center part of the stent, we applied **a hard and long olive tip**, unlike our existing products.
- Stents come in a variety of lengths, giving you a variety of options.



Delivery System 8Fr, 50/180cm  
Order NO. First two\*\* : Diameter(mm) / Last two\*\* : Length(cm)

OTW	TTS	Diameter	Center	A	B	Total
SBP****HCN SBP****HC	SBE****HCN SBE****HC	10, 12	25	15	20	60
				20	25	70
				20	35	80
				20	45	90
				30	45	100
				30	55	110
				40	55	120

