

RadioMate Access Kit - F

(Micro Puncture Set)



Kink-resistant
core-wire available



Smooth transition




Needle bevel indicator

Radiomate Access kits provide the high quality and precision design expected by clinicians worldwide.

Each feature was designed to aid in the access process, allowing for consistent performance every time.


FEATURES AND BENEFITS

- ☑ 4 and 5 French sizes available
- ☑ Variety of wire configurations
- ☑ Standard or Echogenic needles
- ☑ All inclusive packing



ORDERING INFORMATION

Ordering No.	Catheter Sheath Diameter(Fr)	Needle		Guidewire	
		I.D(inch)	O.D(inch)	O.D(inch)	Length(cm)
RMAK0410F	4Fr	0.023	0.032	0.018	40
RMAK0510F	5Fr				



Needle



Gold Tip Wire